# Large Scale Cell Culture Single Pass Tangential Flow Filtration

#### Confidential Client

#### **MANDATE**

Biopharm Engineered Systems LLC (BPES, a CPS Group Company) was retained to Design & Fabricate a Single Pass Tangential Flow Filtration process skid for Confidential Client as part of a large expansion project encompassing around (2) 10,000L bioreactors including both upstream and downstream equipment. Confidential Client specializes in pharmaceutical medicines, vaccines and consumer healthcare products & has the broadest portfolio in the vaccines industry including pediatric, adolescent, adult and traveller vaccines.

#### **PROJECT SCOPE**

BPES's main project scope included merging both a Tangential Flow Filtration (TFF) process and a Single Pass Tangential Flow Filtration (SPTFF) process all onto one skid. Project challenges included interlacing three individual systems onto one skid as well as designing a skid break to accomodate accessibility into the suite. This included meticulous consideration of instrument placement, minimizing tri-clamp placement, designing a hinged ladder & considering routing paths for the breakage of electrical and pneumatic wireways.

## **Process Equipment** (Including but not Limited to)

500L Tulip Vessel

Magnetically Coupled Agitator

2 TFF Feed Pumps

TFF Retentate Vessel CIPR Pump

TFF Recirculation Pump

2-High Ultrafilter Holder (20-40 m2 capacity)

(4) SPTFF holders - (8-20 m2 total capacity)

### Software / Communications

CHARMS I/O





# **Equipment / Services Provided (Including** but not Limited to):

- Mechanical Engineering
- TFF Retentate Vessel & Agitator
- Tangential Flow Filtation Skid
- Single Pass Tangential Flow Filtation
- Filter Housing
- Multi-port block valve engineering
- Hangers and Supports
- Vibrations and Seismic Control
- Hygenic System Specialties
- Piping Systems
- Low Voltage Induction Motors
- Pumps
- **Electrical Engineering**
- Variable Frequency Controllers
- Basic IC Reqs and PCS Controls
- Pneumatic Instrument Tubing
- Instrument Calibration & Commissioning
- Instrument and Control Panels
- Factory Acceptance Testing (FAT)
- Tagging and Identification
- Enhanced Turnover Package
- **GMP** Documentation

# SECTOR OF ACTIVITY

Pharmaceutical and

Biopharmaceutical

Biopharmaceutical

#### CLIENT

Confidential Client

# YEAR OF COMPLETION

2018

## SKID DIMENSIONS

~ 16' x 7' x 12'

## **BPES PROJECT VALUE**

Confidential Value

